The solution for **metal-free** CAD/CAM restorations

MADE IN THE USA

Bicon, LLC
501 Arborway
Boston, MA 02130
USA
tel: 617.524.4443
fax: 617.524.0096
support@bicon.com

Bicon Europe, Ltd.
Unit 4 Ballycummin Village
Ballycummin, Limerick
Ireland
tel: +353 6130 3494
fax: +353 6122-7606
ireland@bicon.com

Subsidiaries\*

Bicon Europe
Limerick, Ireland
ireland@bicon.com
Bicon Italia
Cisterna, Italia
italy@bicon.com

Bicon Vietnam
Ho Chi Minh City, Vietnam
vietnam@bicon.com.vn

\*For complete address information of Bicon subsidiaries and other regional distributors, please visit www.bicon.com

Australian Sponsor

Emergo Australia
201 Sussex Street - Darling Park, Tower II, Level 20
Sydney, NSW 2000 Australia

EN ENGLISH INSTRUCTIONS FOR USE

General

TRINIA™ CAD/CAM discs and blocks are composed of a multi-directional interlacing of fiberglass and resin in several layers. TRINIA™ is milled into the desired shape of a substructure using a dental CAD/CAM system. The restorations are then bonded to tooth structure and/or to implant abutments with adhesive resin cement for use in edentulous or partially edentulous patients.

Indications

TRINIA™ is intended to be used by dental technicians and dentists for making copings, substructures, removable dentures, or frameworks for permanent and transitional anterior or posterior crowns, bridgework, and substructures that can be either cemented or uncemented restorations (e.g. telescopic restorations).

Contraindications

TRINIA™ should not be used in patients with parafunctional habits, such as bruxism.

Warnings and Precautions

- TRINIA™ is supplied in a ready state. This material should not be fired under any circumstances.
• TRINIA™ should not be milled or used at temperatures above 150°C (302°F).
• Do not contaminate TRINIA™ with oils or grease.
• Do not exceed the mechanical tolerances of the device. TRINIA™ has a flexural strength of >300 MPa and a flexural modulus of elasticity of <20 GPa and requires an occlusal thickness of 2.0mm for bars. Failure to observe these parameters may compromise the results achieved with TRINIA™.
• If patients are known to be allergic to any of the ingredients, TRINIA™ restorations should not be used. The processing of TRINIA™ discs and blocks produces dust which may irritate the skin and eyes or cause other health issues. Refer to the Material Safety Data Sheet (available at www.bicon.com).
• Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established, if applicable.

CAUTION: USA Federal law restricts this device to sale by or on the order of a licensed dentist or dental professional.

Preparation

Proper preparation of the tooth/teeth is essential for optimum strength, color, and retention of a milled restoration. TRINIA™ is not water soluble and can be milled with air coolant by following the instructions of the manufacturer of the milling machine. Recommended bur material: Nano-diamond

How Supplied

TRINIA™ is provided in protective packaging with a label containing its dimensions and product code. TRINIA™ is provided in a ready-to-use state. Remove the product from its packaging and check for chips, cracks, or other defects.

Storage and Use

TRINIA™ is designed to be stored in its original packaging and used at room temperature. Do not store the product in heat or in sunlight.

Disposal

The material is non-hazardous and can be disposed of in accordance with regional and local authority requirements.

The solution for metal-free CAD/CAM restorations

Key

Table with 2 columns: Icon and Description. Includes Manufacturer, Date of manufacture, Use by date, EC REP, LOT, REF, Read instructions for use, By prescription only, Caution, Do not use if package is damaged, Keep away from sunlight, Manufacturer, Medical device.

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