





UNIVERSAL ABUTMENT BLANKS



*“Custom abutments with the quality and precision
of the Bicon locking taper connection.”*



Product Guide

Bicon's Universal Abutment Blanks provide for custom abutments while assuring the quality and precision of the Bicon locking taper connection.

- Design custom Bicon abutments with optimal gingival contours
- Use your own in-house milling machine
- Increase your efficiency
- Assure maximum titanium contact with gingival sulcus, especially for less than ideally positioned implants

Material: Ti6Al4V-ELI

Connection: Bicon Locking Taper




2.0mm Post Universal Abutment Blanks

Part No.	Diameter	Series
260-741-700	10mm	UA4
260-744-700	14.2mm	UA4



2.5mm Post Universal Abutment Blanks

Part No.	Diameter	Series
260-541-700	10mm	UA4
260-544-700	14.2mm	UA4
260-551-700	10mm	UA5
260-554-700	14.2mm	UA5



3.0mm Post Universal Abutment Blanks

Part No.	Diameter	Series
260-651-700	10mm	UA5
260-654-700	14.2mm	UA5
260-661-700	10mm	UA6
260-664-700	14.2mm	UA6

