






PILOT DRILL STOP KIT



"Allows experienced and novice clinicians alike to confidently drill a 2.0mm pilot osteotomy to a precise depth."



Highlights


When the anatomy permits, Bicon's 2.0mm diameter Pilot Drill Stop Kit allows experienced and novice clinicians alike to confidently drill a 2.0mm pilot osteotomy to a precise depth.

- Seven stops ranging from 5.0mm–11.0mm
- Self-contained in compact surgical tray
- Stainless steel with titanium-nitride coated tip

Material: Stainless Steel

Drills

Part No.	Diameter	Stop
101-105	2.0mm	5.0mm
101-106	2.0mm	6.0mm
101-107	2.0mm	7.0mm
101-108	2.0mm	8.0mm
101-109	2.0mm	9.0mm
101-110	2.0mm	10.0mm
101-111	2.0mm	11.0mm



Kits

Part No.	Description
101-099	2.0mm Pilot Drill Stop Kit
104-999	2.0mm Pilot Drill Stop Tray



2.0mm Pilot Drill Stop Kit